

Attaching Fluorescein to a single cysteine in proteins

Labeling buffer (45 ml):

10 mM K-Hepes (pH 7.0)	1M stock	0.45 mL
100 mM KCl	3M stock	1.5 mL
20 μ M GDP*	100 mM stock	9.0 μ L
Water		43.0 mL

Note: GDP is needed only for EF-Tu or RF3

5 mM Iodoacetamidofluorescein (FW =515):

Freshly dissolve 1.03 mg in 400 μ L of Dimethylformamide (DMF). Do not store, discard unused reagent.

250 μ M β -Mercaptoethanol:

Mix 2 μ l (14.4 M) mercaptoethanol stock with 498 μ l of water = 57.6 mM stock.

Mix 2 μ l (57.6 mM stock) with 458 μ l of water = 250 μ M stock.

Protein is stored in Labeling Buffer: 10 mM HEPES-KOH (pH 7.0), 100mM KCl, 20 μ M GDP*.

Labeling Reaction:

125 μ M Protein	100 μ M final conc.	160 μ l
250 μ M β -mercaptoethanol	25 μ M final conc.	20 μ l
5 mM Iodoacetamidofluorescein	500 μ M final conc.	20 μ l

Add Labeling buffer to 200 μ l final volume.

1. Mix EF-G proteins (100 μ M) and limiting β -mercaptoethanol (β -ME; 25 μ M) first and incubate at room temperature for 1 hour to reduce the cysteines. Then add the Iodoacetamidofluorescein to start the labeling reaction.
2. Incubate reaction in the dark (20°C, overnight).
3. Stop reaction by adding 1.4 μ l of 14.4 M β -ME (final conc. is 100 mM).
4. Purify EF-G by passing through a Sephadex G-25 or BioRad P2 column (make sure that the molecular weight cut-off is appropriate for the column).
5. Analyze a small amount by SDS-PAGE. The labeled protein will fluoresce under UV light.

6. Dialyze against storage buffer, and store at -80°C in small aliquots.

The ratio R of FL attached to EF-G varied between 0.6 and 0.8, as determined by the protein's absorbance of light at 280 and 492 nm, and the equation:

$$R = (A_{492} \times \epsilon_{EF-G}) / (A_{280} \times \epsilon_{FL}),$$

where $\epsilon_{FL} = 78,000 \text{ M}^{-1}\text{cm}^{-1}$ and $\epsilon_{EF-G} = 64,282 \text{ M}^{-1}\text{cm}^{-1}$.

NOTES:

-The dye (FL) is not stable after dissolving in DMF. So dissolve only small quantities by weighing few specks on weighing paper and adding DMF to appropriate volume.

-It may be necessary to purify the protein by gel filtration after labeling.