

#### **Protocol 4: Deprotection and Purification of Synthetic DNA**

*CAUTION: Wear gloves, use RNAase free materials.*

1. Carefully open the column and transfer CPG beads into a 1.5 ml screw-cap tube.
2. Add 500  $\mu$ l of Conc. Ammonium hydroxide. Close the tubes tightly and seal with parafilm.
3. Incubate at 55 °C overnight. Cool to -20 °C for 10 minutes. Spin at max for 1 min and transfer supernatant to new 1.5 ml eppendorf tubes.
4. Dry above mix to a white powder in Speed Vac with Heat. Takes 4-5 hours.
5. Resuspend dried pellet in 100  $\mu$ l TE, pH 7.5.
6. Further purify the DNA on a denaturing polyacrylamide gel. Load about 30  $\mu$ l DNA + 30  $\mu$ l 2XSLB mix/lane. Follow standard procedures for DNA identification by UV-shadowing and passive elution from gel slice.