

Protocol: Double Stranded DNA Sequencing

I. Preparing Template:

1. Plasmid DNA = 10 μ L (3-5 ug)
(2M) NaOH = 2 μ L
Total vol = 12 μ L
2. Incubate at 37 C for 20 min
3. Add 10 μ L of (3M) NaAc (pH5.2) + 78 μ L of TE
4. Extract with 100 μ L of phenol/chloroform mix (pH 7.0. BRL)
5. Extract with 100 μ L of chloroform
6. Add 300 ul of ice-cold 100% ethanol. Precipitate at -80 C for 10 min.
7. Spin in microfuge for 30min at 4 C
8. Wash with 100 ul of 70% ethanol. Dry pellet. Resuspend in 6.0 μ l of water.

II. Annealing Reaction:

1. Add 2.0 μ l of Sequenase reaction buffer and 10 picomoles (1 μ l) of primer to above.
2. Incubate at 70 C for 2 min. Cool slowly to room temp. Place on ice.

III. Labeling Reaction:

1. Dilute the 5X Labeling mix with water (example: 4 ul of LM + 16 ul of water). Dilute Sequenase enzyme 1:8 with Enzyme Dilution Buffer (example 2.5 l of enzyme + 17.5 ul of enzyme dilution buffer). Leave both on ice. Use immediately.
2. Extension reaction:

Template+Primer mix	= 9.0 μ l
100 mM DTT	= 1.0 μ l
Diluted Labeling mix	= 2.0 μ l
[α -32P] dTTP	= 1.0 μ l
Diluted Sequenase	= 2.0 μ l
Total vol	= 15.0 μ l
3. Mix and incubate at room temp for 5 min.

IV. Termination Reaction:

1. Label four tubes G, A, T, C. Add 2.5 μ l of ddGTP, ddATP, ddTTP and ddCTP termination mix to the four tubes respectively. Incubate at 37 C for 2 min.
2. Add to each tube 3.5 μ l of the Labeling Reaction mix.
3. Incubate at 37 C for 5 min.
4. Add 4 μ l of 2X Sequencing Loading Buffer (2XSLB) to stop the reaction.
5. Store on ice until ready to load on 6% sequencing gel.

6. Heat samples at 90 C/2min. Place on ice. Load immediately 1.5 μ l/ lane.

NOTES:

To sequence close to the primer dilute the 5X Labeling mix to 0.25X (20-fold) and add 2.0 μ l to the Labeling Reaction. Also add 1 μ L of supplied Mn Buffer to the Labeling reaction.