

Determining the efficiency of aminoacylation reaction

1. Properly fold the tRNA in 30 mM HEPES (pH 7.0) and 30 mM KCl by heating to 70 °C for 5 min after which 15 mM MgCl₂ is added and the reaction is then allowed to cool to room temperature over a period of 30 min by transferring to metal tray containing hot water.
2. Perform small scale aminoacylation reaction with [³⁵S] methionine or [¹⁴C] phenylalanine. Example for *E. coli* tRNA^{fmet} is described below:

30 mM Hepes (pH 7.0)

30 mM KCl

15 mM MgCl₂

8 mM ATP (pH 7.0)

10 mM DTT

20 μM Methionine

0.05 μM [³⁵S] methionine (specific activity = 5,000 cpm/picomoles of methionine)

5 μl of MetRS or S-100 extract

10 μM *E. coli* tRNA^{fmet}

Final volume = 100 μl (make enough for duplicates and control reactions)

3. Incubate at 37 °C for 20 min.
4. For calculating specific activity of the reaction mix: Directly spot 2 μl of the reaction mix from step 2 onto filter and count without washing.
5. Take 50 μl aliquot of the aminoacylation reaction and mix with 50 μl of carrier (1 μg/ml BSA, 1M sodium acetate pH 5.2) and 1 ml of precipitation mix (2.5% trichloroacetic acid in 50% ethanol).
6. Place on ice for 20 minutes.
7. Filter the above mix through Whatman GF/C filters using the Millipore Filtration setup.
8. Wash the filter with 5 ml of 5% trichloroacetic acid. REPEAT this step 4 times.
9. Wash the filter with 5 ml of 95% ethanol. REPEAT this step 4 times.
10. Dry the filter under heat lamp for 10 minutes.
11. Add 4 ml of scintillation cocktail and count in a scintillation counter for the correct isotope.

Controls:

Perform reaction without any tRNA for calculating background.