

### Protocol 6: Fe-BABE Modification of RNA with 5'-Phosphorothioate.

1. Mix 20  $\mu\text{l}$  of 28 mM freshly made ammonium iron (II) sulfate hexahydrate (Aldrich#20350-5) with 25  $\mu\text{l}$  of 28 mM BABE.
2. Adjust pH to 3-4 using 1-2  $\mu\text{l}$  diisopropylethylamine (pH 9). Use pH paper to verify pH.
3. Incubate at room temperature for 1 hour.
4. Combine:

5'-GMPS-RNA	=	$\mu\text{l}$	(800 picomoles)
12 mM Fe(II)-BABE mix (from step 3)	=	5.0 $\mu\text{l}$	
400 mM Potassium phosphate (pH 8.5)	=	1.0 $\mu\text{l}$	
Water	=	$\mu\text{l}$	
		-----	
		10.0 $\mu\text{l}$	Final volume
		-----	
5. Incubate at 37  $^{\circ}\text{C}$  for 1 hour.
6. Add 1.0  $\mu\text{l}$  of 50 mM EDTA. Incubate for 10 minutes at 37  $^{\circ}\text{C}$ .
7. Add 10  $\mu\text{l}$  of 5M ammonium acetate, 1  $\mu\text{l}$  of 10mg/ml glycogen and 80  $\mu\text{l}$  of water.
8. Extract with 100  $\mu\text{l}$  of water-saturated Phenol. Transfer supernatant to new tube.
9. Repeat Phenol and then chloroform extractions. Transfer supernatant to new tube.
10. Add 300  $\mu\text{l}$  of 100% ethanol. Place in  $-80^{\circ}\text{C}$  for 1 hour or  $-20^{\circ}\text{C}$  overnight.
11. Spin tubes in microfuge at max speed for 30 minutes at 4  $^{\circ}\text{C}$ .
12. Carefully remove ethanol without disturbing the pellet. Dry the pellet in SpeedVac.
13. Resuspend the RNA pellet in water at 25 picomoles/  $\mu\text{l}$  concentration.
14. Store Fe-BABE-RNA at  $-20^{\circ}\text{C}$  and use it within 2-3 weeks for tethered hydroxyl radical probing experiments.

**Note: We currently have following stocks of BABE:**

- (a) 28 mM BABE in 20  $\mu\text{l}$  aliquotes synthesized on 8/7/95
- (b) 49.6 mM BABE in 20  $\mu\text{l}$  aliquotes synthesized on 4/22/97.