

Ribosome Filter Binding Assay (Xinying Shi)

- **tRNA Dephosphorylation:**

(100pmol/ul) tRNA	25ul	(2500pmol)
10x buffer 3(NEB)	50ul	
(1u/ul) CIP AK Phosphatase (Amersham)	25ul	(1u/100pmol)
<u>water</u>	<u>400ul</u>	
	500ul	x2

1. 37C 60min,
2. Add 50ul 3M NaOAc,pH 5.2
3. 1x phenol extraction
4. 2x chloroform extraction
5. 2.5x EtOH precipitate
6. 4C 30min,
7. 16000g 30min
8. Resuspend in 40ul water

- **tRNA Labeling**

(100pmol/ul) Dephosphorylated tRNA ^{Phe}	5.0ul
³² P- γ -ATP	5.0ul
10 x PNK	2.5ul
(10u/ul) PNK(NEB)	2.0ul
<u>water</u>	<u>10.5ul</u>
	25.0ul

1. 37C 30min,
2. Add 25ul 2x sequence loading buffer
3. Heat at 90C for 2min,then place on ice for 10min
4. Purify on 10% dPAGE gel, 30W 2h 30min
5. Exposure on film for 5min and develop for 2min (place four fluorescence marker on each conner of the gel)
6. Put the film on the gel and cut out the band
7. Add 300ul RNA elution buffer and shake O/N
8. Transfer buffer to a new tube and wash gel with 300ul elution buffer.
9. 2x chloroform extraction (something white in the interface)
10. Add 2x 750ul EtOH, -20C o/n or 1h
11. 16000g, 30min
12. Resuspend in 30ul water.
13. Add cold tRNA^{Phe}
14. Measure A₂₆₀ of hot&cold tRNA
15. Measure radioactive intensity by scintillation counting.
16. Adjust to 1pmol= 1000cpm by adding cold tRNA.

Note: Store labeled tRNA at -20C for at most 1 week! Use fresh ATP!

Buffer: C₈₀M₁₀N₁₅₀ or C₈₀M₂₀N₁₅₀

	Stock solution	1x	5x
80mM K-cacodylate	1M,pH7.5	4.0ml	400ul
10mM MgCl ₂	1M	0.5ml	50ul
150mM NH ₄ Cl	5M	1.5ml	150ul
water		44.0ml	400ul
		50.0ml	1000ul

- Prepare labeled tRNA.
- Prepare 50ml 1x CMN buffer and store at 4C o/n
- Clean and set up filtration apparatus(Millipore), adjust flow rate (drip, not in a line)
- Refold tRNA in 10mM MgCl₂ at 60C for 3min,then slowly cool down to RT and stay at RT for 5min, then add 5x CMN and dilute to 1x
- Soak nitrocellulose filter(2.5cm,0.45uM, HWAP, Millipore, cat#) in 1x CMN for at least 1h at 4C

Note: Don't use buffer containing polyamine. It will increase the background.

Reaction condition for P site binding:

30s or 70s	10pmol	10ul	reactive at 42C,10min---37C---37C,10min
mRNA	20pmol	5ul	
tRNA ^{Phe}	40pmol	10ul	refolded
	total volume =25ul		

Do duplicate for each reaction.

Need negative control (no ribosome)

1. Refold tRNA at 60C for 3min, slowly cool to RT, and stay for 5min.
2. Reactive at 42C 10min, slowly cool to 37C, stay at 37C for 10min.
3. Associate ribosome, tRNA and mRNA at 37C, 30min.
4. Place on ice.
5. Wash nitrocellulose filter with 1ml cold CMN buffer.
6. Dilute 25ul sample to 1ml with 1x CMN and add to the filter immediately.
7. Wash filter with 2ml buffer.
8. Place filter in order and dry under the heating lamp for 3 min.
9. For counts of each pmol, add 2ul of tRNA^{Phe} to Whatman filter circle (2.3mm, grade 3)
10. Place filter into a marked scintillation vial.
11. Add 4ml scintillation solution (Fisher).
12. Determine radioactivity by scintillation counting.
13. Calculate the binding efficiency:

$$\frac{\text{counting} - \text{background}}{\text{counts per pmol}}$$

Note: About 10% background.