

Protocol for Assigning R_P vs S_P to Separated Phosphorothioate Diastereomers

Ref: Slim and Gait, Nucleic Acids Res 1991, **19** pp1183-1188

1) Set up the reaction:

100 pm/μl RNA	5.90	100 pm/μl RNA	5.90
100 mM Tris-Cl pH 8.5	5.00	1 M Tris-Cl pH 8.5	5.00
1 μg/μl P1 Nuclease (~5U)	5.00	3 mM DTT	5.00
2 μg/μl Alkaline phosphatase (~20U)	5.00	3 mM MgCl ₂	5.00
MilliQ H ₂ O	29.10	5 μg/μl SVPD (~0.5U)	5.00
		2 μg/μl Alkaline phosphatase (~20U)	5.00
		MilliQ H ₂ O	19.10

Total volume = 50 μl

2) Incubate 12-16 hrs at 37 °C

3) Spin down condensate and incubate at 90 ° C for 3 minutes

4) Place on ice or store at -80 ° C until ready to run on HPLC

Buffer A: 20 mM triethylammonium acetate (TEAA) pH 7

Buffer B: 60 % Buffer A, 40% acetonitrile

Method: reverse phase, using Hamilton PRP-1, 7 mm * 305 mm column

HPLC settings:

- 1) Use 100 μL injection loop, ~50 μL injection volume
- 2) Flow rate of 1.5 mL/min (pressure is around 1.000)
- 3) Chart rate of 30 cm/hr
- 4) Gradient of 5% B → 65% B, 30 minutes
- 5) Set chart recorder to 2 mV on most sensitive channel

Materials:

2 M Triethylammonium acetate pH 7: Fluka cat # 09748

Acetonitrile: Fisher Optima Grade, cat # A996

Nuclease P1: Roche cat # 236 225, 1 mg dissolved in 1 mL 30 mM Na-acetate pH 5.3 to about 1 U/μl and stored at 4 ° C

Snake Venom Phosphodiesterase I: Worthington Biochem cat # LS 003926, 1 mg dissolved in 100 mM Tris pH 8.5 to about 0.1 U/μl, and stored at -20 ° C

Alkaline Phosphatase: Roche cat # 405612, 2 mg dissolved in 1 mL AP Buffer to 4 U/μl, and stored at 4 ° C

AP Buffer: Volume = 1mL

4 M NaCl	25.00
1 M MgCl ₂	10.00
100 mM ZnCl ₂	10.00
1 M Triethanolamine pH 7.6	30.00
MilliQ H ₂ O	925.00

Protocol written by: Steve Phelps, 2004