

## Protocol: Preparation of E. coli S-100 Extract for Luciferase Assay (Jill Thompson)

**IMPORTANT:** Everything should be done at 4 °C! Use ice-cold Buffers, Rotors, etc. Use RNase free lab wares and reagents.

### Materials

#### Wash Buffer

Stock	Final Conc	Volume (1L)
1M HEPES-KOH, pH 7.5	10 mM	10 ml
1M MgCl <sub>2</sub>	10 mM	10 ml
2.5M NH <sub>4</sub> Cl	1000 mM	400 ml
Mercaptoethanol	5 mM	350 ul
Water		580 ml

#### S-30 Buffer

Stock	Final Conc	Volume (50 ml)
1M HEPES-KOH, pH 7.5	10 mM	500 ul
1M MgCl <sub>2</sub>	10 mM	500 ul
2.5M NH <sub>4</sub> Cl	50 mM	1000 ul
Mercaptoethanol	5 mM	18 ul
Water		48 ml

#### S-30 Buffer + 10% glycerol

Stock	Final Conc	Volume (1L)
1M HEPES-KOH, pH 7.5	10 mM	10 ml
1M MgCl <sub>2</sub>	10 mM	10 ml
2.5M NH <sub>4</sub> Cl	50 mM	20 ml
100% Glycerol	10%	100 ml
Mercaptoethanol	5 mM	350 ul
Water		860 ml

#### 100 mM Glycine-KOH, pH 9.2:

75 mg of glycine (F.W = 75) in 10 ml water. First pH with KOH to 9.2 before final volume is made to 10 ml.

#### 25 mM EGTA, pH 7.0:

95 mg of EGTA (F.W = 380.4) in 10 ml water. pH with KOH.

#### 50 mM CaCl<sub>2</sub>:

74 mg CaCl<sub>2</sub>·2H<sub>2</sub>O in 10 ml water

**Procedure:**

1. Inoculate FOUR 5 ml LB with MRE600 glycerol stock. Grow overnight at 37 °C in shaker.
  2. Inoculate TWO 500 ml LB with 5 ml of saturated, overnight MRE600 culture. Grow at 37 °C in a shaker until A650 is 0.5. (Takes approximately 2-3 hours)
  3. Transfer to 500 ml centrifuge bottles.
  4. Spin at 8000 RPM for 15 min at 4 °C in BECKMAN JA10 rotor.
  5. Resuspend each cell pellet in 15 ml of Wash Buffer. Spin as in step 4. Resuspend again in 15 ml of S-30 Buffer.
  6. Pool all the resuspended cells in one centrifuge bottle. Spin at 8000 RPM for 15 min at 4 °C in BECKMAN JA10 rotor. You can store the cell pellets at this stage at -80 °C.
  7. Weigh the cell pellet. Resuspend the cell pellet in S-30 Buffer + 10% glycerol on ice:  
Add 1.5 ml buffer per Gram of cells
  8. Lyse the cells using a French Press.at 10,000 PSI into a red-capped ultra centrifuge bottle (BECKMAN Cat# 355618) on ice. (Note: Before use, rinse the French press cylinder with water, then with S-30 buffer and cool the cylinder to 4 C)
  9. Centrifuge at 20,000 RPM for 30 min at 4 C in Ti70 rotor.
  10. Transfer the supernatant to a new ultra centrifuge bottle (discard the pellet). Centrifuge again at 20,000 RPM for 30 min at 4 C in Ti70 rotor.
  11. Pool the supernatant and measure the A260, which should be about 150-250 units/ml. (Save about 20 ul at -80 C for assay).
  12. Nuclease Treatment:  
Dissolve 15,000 units of Micrococcal nuclease (Boehringer cat# 107921) in 200 ul of buffer (50 mM glycine-KOH, pH 9.2; 5 mM CaCl<sub>2</sub>).
- To 1 ml of S-30 supernatant from step 11, add 20.5 ul of 50 mM CaCl<sub>2</sub> (final conc 1 mM) and 10 ul of Micrococcal nuclease (750 units). Incubate at 30 C for 20 min. (Note: I recovered 5.5 mls of S-30, so added 113 ul of 50 mM CaCl<sub>2</sub> and 55 ul of nuclease).
13. Inactivate the nuclease by adding 86.65 ul of 25 mM EGTA-KOH, pH 7.0 (final conc 2 mM). (save 20 ul for assay). (I added 440 ul of 25 mM EGTA-KOH).

14. Centrifuge at 50,000 RPM for 2 hours at 4 C in Ti70 rotor. Note: The tube will collapse at 50K unless it contains 16 ml minimum volume.
15. Remove the supernatant. Dialyze against S-30 Buffer + 10% glycerol. Use 6-8K cut-off dialysis tubing. Dialyze TWO times against 100 volume for 45 minutes each at 4 C.
16. Make 50 ul aliquots and store at -80C. (About 100 aliquotes)