

Protocol for Separating Phosphorothioate Diastereoisomers via HPLC

Method 1: Reverse phase, using Waters C18 μ -Bondapak, 3.9 mm * 300 mm column
Ref. Sam Butcher communication; Slim and Gait, Nucleic Acids Res 19 pp1183-1188, 1991

Buffer A: 100 mM ammonium acetate, pH 6.8

Buffer B: 20% A, 80% acetonitrile

HPLC settings:

- 1) Use 100 μ L injection loop, ~50 μ L injection volume
- 2) Run about $\frac{1}{4}$ of the deprotected RNA prep from a 0.4 μ m synthesis (~150 μ g)
- 3) Flow rate of 1.5 mL/min (pressure is around 1.000)
- 4) 0-8 min, 0% B
8-11 min, 0-5% B
11-41 min, 5-7% B (optimized for 19mer RNA)
41-49 min, 20% B
- 5) Chart rate of 30 cm/hr, 20 mV on the medium sensitivity channel
- 6) Run at room temperature

Method 2: Anion exchange, using Dionex DNAPac PA-100, 9 mm * 250 mm column
Ref. Dionex user's manual

Buffer A: 100 mM ammonium acetate, pH 8; 370 mM KCl; and 2% acetonitrile

Buffer B: 100 mM ammonium acetate, pH 8; 700 mM KCl; and 2% acetonitrile

HPLC settings:

- 1) Use 100 μ L injection loop, ~50 μ L injection volume
- 2) Run about $\frac{1}{4}$ of the deprotected RNA prep from a 0.4 μ m synthesis (~150 μ g)
- 3) Flow rate of 1.5 mL/min (pressure is around 0.500)
- 4) Gradient of 0 \rightarrow 56% B over 30 min (optimized for 19mer RNA)
- 5) Chart rate of 30 cm/hr, 20 mV on the medium sensitivity channel
- 6) Heat column to 65 $^{\circ}$ C, then equilibrate in buffer A

Materials:

Ammonium acetate: Fisher, HPLC grade, cat. # A639

KCl: Sigma, 99% pure, cat. # P-9541

Acetonitrile: Fisher, Optima grade, cat # A996

Notes:

- 1) Scale down to less than 100 μ g if resolution is poor
- 2) Make a more shallow gradient around the %B that you know the oligo will elute
- 3) Always filter samples through a 0.22 μ m filter before injection
- 4) Use Waters "SepPak" Cartridge for desalting, but beware of high acetonitrile content in your eluate after reverse phase HPLC
- 5) If Method 1 doesn't give good separation, try Method 2

Protocol written by: Steve Phelps, 2004