

Purification of U1A Tagged 50S subunits (Jason Feinberg) version 2

This protocol is a morph of techniques from protocols provided by the Strobel lab and the Noller lab. It has been optimized for time and yield. Variations include usage of French press, deletion of pelleting spins and deletion of zonal centrifugation. It is important to note that purity of His-tagged 50S depends greatly on removing the 30S as can be seen in poison primer experiments. (J.S.F 2006)

- Inoculate 250 ml LB containing chloramphenicol and ampicillin with scrape from His₆-U1A cell stock.
- Grow/shake o/n @ 30°C
- Add 5 ml of o/n culture to (one or more) 1 l LB containing chloramphenicol and ampicillin **pre-warmed to 42°C**.
- Grow/shake at 42°C to a final OD₆₀₀ of 0.5
- Ice for 20 min.

- Pellet cells at 6-8k rpm in Beckman JA-10.
- Wash (resuspend and spin as above) and resuspend cells in 15 ml of Cell Lysis Buffer.
- Lyse cells in French press (2 passes ~18000 psi)
- Add 10 µl RQ1-DNAse
- Pellet cell debris at 15k rpm in Beckman JA-17
- Transfer supernatant to fresh tubes and spin again @ 15k rpm

- Using “red capped” Ti-70 centrifuge tubes put 9 ml sucrose cushion under 14 ml of supernatant from previous spin (adjust volume w/wash buffer of supernatant for even # tubes).
- Spin C

- Wash pellets 1 ml of Washing Buffer and resuspend in 10 ml Washing Buffer.
- Pool resuspended pellets and place into “red capped” centrifuge tubes with a 1 ml sucrose cushion. Spin D.

- Wash pellets with 1 ml of **LOW Mg⁺²** buffer and resuspend pellets in 10 ml of **LOW Mg⁺²** buffer. (you now have a crude ribosome slurry with dissociated 30S and 50S)

Take OD₂₆₀ of ribosomes and calculate nmoles.

$$((\text{OD}_{260} \times \text{dil. Factor}) / 39.15) \times \text{volume in mls} = \text{nmoles of ribosomes}$$

For every 40 nmoles of ribosomes, use 1 ml Ni-NTA beads

- Prepare Ni-NTA beads by rinsing with 2x in MiliQ H₂O in a 50 ml falcon. Then rinse in 1x PD buffer.
- Add ribosome slurry to Ni-NTA beads using the following formula:

For every ml of beads there should be 7.5x volume of solution. So: 5 ml of beads have a final volume of beads + ribosomes + buffer = 35mls

So:

10 mls ribosomes in LOW Mg²⁺ buffer
 5 mls beads in 1X PD buffer
20 mls 1X PD buffer
 35 mls final volume

- Mix ribosomes with beads gently rocking for 20 min @ 4°C
- Pour over a pre-chilled column
- Collect flow through and re-pour over column
- Wash with at least 20 x column volumes of 1x PD buffer
- Drain until just above beads and add 20ml of 1xEB
- Collect eluate and put into ultra “red capped” tubes. Balance with Ribosome Storage Buffer.
- Spin 3.5 hours at 60K
- Resuspend in ~0.1-1 ml Buffer A depending on pellet size
- Take O.D. and calculate concentration of 50S
- Aliquot into eppendorf tubes, flash freeze with liquid nitrogen and store at -80°C

<u>Cell Lysis Buffer</u>	<u>Stock</u>	<u>500 ml</u>
20 mM Tris-HCl (pH 7.6)	1M	10 ml
100 mM NH ₄ Cl	5M	10 ml
10.5 mM Mg acetate	1M	5.25 ml
0.5 mM EDTA	0.5M	0.5 ml
3 mM β-mercaptoethanol	14.4	10.5 μl
<u>Sucrose Cushion</u>	<u>Stock</u>	<u>200 ml</u>
20 mM Tris-HCl (pH 7.6)	1M	4 ml
100 mM NH ₄ Cl	5M	20 ml
10.5 mM Mg acetate	1M	2.1 ml
0.5 mM EDTA	0.5M	0.2 ml
Sucrose		75.24 g
3 mM β-mercaptoethanol	14.4	4.2 μl
<u>Washing Buffer</u>	<u>Stock</u>	<u>500 ml</u>
20 mM Tris-HCl (pH 7.6)	1M	10 ml
500 mM NH ₄ Cl	5M	50 ml
10.5 mM Mg acetate	1M	5.25 ml
0.5 mM EDTA	0.5M	0.5 ml
7 mM β-mercaptoethanol	14.4	10.5 μl
<u>Low Mg⁺² Washing Buffer</u>	<u>Stock</u>	<u>500 ml</u>
20 mM Tris-HCl (pH 7.6)	1M	10 ml
500 mM NH ₄ Cl	5M	50 ml
1.5 mM Mg acetate	1M	0.74 ml

0.5 mM EDTA	0.5M	0.5 ml
7 mM β -mercaptoethanol	14.4	10.5 μ l
<u>Pull Down Buffer</u>	<u>Stock</u>	<u>200 ml</u>
50 mM Tris-HCl (pH 7.6)	1M	10 ml
5 mM Imidazole	1M	1ml
70 mM NH_4Cl	5M	2.8 ml
1 mM MgCl_2	1M	200 μ l
10 μ g/ml total tRNA	10 mg/ml	200 μ l
10 μ g/ml Heparin	100 mg/ml	20 μ l
<u>Elution Buffer</u>	<u>Stock</u>	<u>20 ml</u>
50 mM Tris-HCl (pH 7.6)	1 M	1 ml
70 mM NH_4Cl	5 M	0.7 ml
30 mM KCl	3 M	0.2 ml
10 mM MgCl_2	1 M	2 ml
200 mM Imidazole	1 M	4 ml
<u>Ribosome Storage Buffer</u>	<u>Stock</u>	<u>100 ml</u>
50 mM Tris-HCl (pH 7.6)	1 M	5 ml
70 mM NH_4Cl	5 M	3.5 ml
30 mM KCl	3 M	1.5 ml
7 mM MgCl_2	1 M	0.7 ml
DTT	1 M	100 μ l
0.5 mM EDTA	0.5 M	100 μ l

Spin	mls Sucrose	Spin Option
A	x	JA-20 rotor 10 min at 3000 rpm Then 30 min at 9000 rpm
B (4 tubes)	x	Ti-70 30 min at 16K
C (8 tubes)	9 ml	8.6h @ 50K 13.5h @ 40K 15h @ 38K 16.5h @ 36K 17.5h @ 35K
D (4-8 tubes)	1 ml	4.4h @ 65K 4.8h @ 62K 5.2h @ 60K 5.5h @ 58K 6.2 @ 55K 6.7h @ 53K 7.4h @ 50K 8h @ 48K 9h @ 45K 15h @ 35K