

Protocol 2: In vitro Transcription of tRNA

Step1: Digest the plasmid encoding tRNA with BstNI

Water	xx μ l
10X NEB2 buffer	10 μ l
1 μ g/ μ l BSA	10 μ l (not 100X)
BstNI (NEB)	10 μ l
Plasmid	xx μ l (80 μ g)
	<hr/>
	100 μ l
	<hr/>

1. Incubate at 60 °C for 2 hours.
2. Add 10 μ l of 0.5 M EDTA, pH8.0 and 10 μ l of 3M Sodium acetate pH 5.2.
3. Add 120 μ l of phenol-chloroform-isoamylalcohol mix, extract by vortexing rapidly for 30 sec. Spin in microfuge for 1 min at max speed.
4. Transfer upper aqueous phase to new tube. Add chloroform and repeat step 3.
6. To aqueous phase add 350 μ l of 100% Ethanol. Precipate DNA at -80 °C for atleast 30 min or -20 °C for 12 hours.
7. Spin tubes at 4 °C for 30 minutes.
8. Carefully remove and discard ethanol without disturbing the DNA pellet.
9. Add 500 μ l of 70% ethanol, spin at max for 30 seconds. Carefully remove all the ethanol. Dry pellet for 5 minutes in the SpeedVac without heat.
10. Resuspend in 50 μ l of TE.

Step2: Reaction mix for 500 μ l transcription.

Water	223 μ l	
10X Trx Buffer (New)	50 μ l	
50 mM GMP	100 μ l	(10 mM Final Conc.)
50 mM DTT	50 μ l	
25 mM AUGC (NTP mix)	40 μ l	(2 mM Each Final Conc.)
T7 RNA Polymerase (JF prep)	10 μ l	
Digested Plasmid	25 μ l	(40 μ gs)
$[\alpha\text{-}^{32}\text{P}]$ UTP (Optional)	2 μ l	
	500 μ l	

Follow Protocol 1: Steps 1-21.