

Storage Phosphor Screens

Screens for storage phosphor imaging systems

Storage Phosphor Screens are used in conjunction with a storage phosphor system, such as Typhoon™ Variable Mode Imagers, Storm™ Gel and Blot Imaging Systems, and PhosphorImager™ instruments. Storage Phosphor Screens retain energy from beta particles, X-rays, and gamma rays, and require exposure times that are 50–90% shorter than that of conventional film. Upon laser-induced stimulation, light is emitted from the storage phosphor screen in proportion to the amount of radioactivity in the sample.

The resulting digital image allows for quantitation of subtle signal intensity differences over the wide dynamic range of the system and for production of publication-quality images. Latent images are erased with the ImageEraser, a white light source. Screens are reusable and are not degraded by repeated exposure to laboratory levels of radioactivity.

Screen types

General purpose

General purpose screens are reliable for a wide variety of applications and can be used with ^{32}P , ^{125}I , ^{35}S , ^{33}P and ^{14}C . The durable cellulose acetate coating and the phosphor layer formulation makes the GP screen ideal for ^{32}P and ^{125}I detection and quantitation. This is the screen of choice for ^{32}P Northern blots, Southern blots, ^{125}I Western blots, and gels.

Tritium

To detect the weak energy of the ^3H signal, Tritium (TR) screens are constructed without the protective cellulose acetate overlay. The maximum energy range of ^3H is 5 microns and that signal can be further attenuated by a thick sample or sample media.

A sample prepared with a fluorographic reagent is not recommended. Best results are obtained when the ^3H signal is on the surface of the sample and available to penetrate the screen (for example, tissue sections on glass and whole body autoradiography).

Screen styles

Two types of screen configurations are available. Mounted screens are permanently fixed to an aluminum backing plate, and they require a specially designed exposure cassette. Unmounted screens have no additional backing and can be used with a standard exposure cassette.

PhosphorImager Systems (other than Storm and Typhoon) require the use of mounted screens. Storm and Typhoon Systems, with their lid pressing mechanism, can use mounted or unmounted screens.

Detection threshold

The lower limit of detection for a 1-h exposure is less than 2 dpm/mm² for ^{14}C using the GP screen scanned at 200 μm.

Excitation

Creation of excited (stored) electrons occurs when the screen is exposed to wavelengths shorter than 380 nm or to α , β , or γ radiation.

Erasure

Approximately 80% of stored information is released upon scanning with the Phosphor Imager. The remaining signal must be erased before the Storage Phosphor Screen is reused. As shown by the stimulation spectrum (Fig 3), the most efficient wavelengths of erasing light are between 475 and 650 nm.

For example, using a 500 W photoflood tungsten light bulb and a yellow filter (Wratten 16), an exposure of 10 J/cm² will leave a residual signal level of less than 10⁻⁵ of the original exposure level.

Screen selection guide

Instrument compatibility		PhosphorImager	Storm and Typhoon Systems
		Models: 400, 425, 445, PSI and PSF	Models: 820, 830, 840, 860, 8600, 92X0, and 94X0
TYPE		GP (mounted)	GP (unmounted)
Small Phosphor Screen 20 x 25 cm	screen & cassette	63-0034-89	63-0034-86
	screen only	63-0034-88	63-0034-87
Large Phosphor Screen 35 x 43 cm	screen & cassette	63-0034-82	63-0034-79
	screen only	63-0034-81	63-0034-80
Tritium Screen 19 x 24 cm	screen only	63-0035-50	63-0035-49
	screen kit	63-0035-51 (1 mounted, 4 unmounted)	

GP (general purpose screen) – for use with ¹²⁵I, ³²P, ³³P, ³⁵S and ¹⁴C

Performance properties of storage phosphor screens.




Fig 1. Screen resolution.




Fig 2. Latent image retention. Greater than 90% of the stored energy is available two minutes after exposure, and greater than 50% remains up to 24 h.




Fig 3. Stimulation spectrum. Stimulation is the process of releasing the stored information. Phosphor Imaging Systems use a red laser light source to deliver greater than 2 mW to the sample.




Fig 6. Ambient radiation background effect. Screens were erased and then stored at room temperature in the exposure cassette prior to scanning. Ambient radiation is assumed to be environmental gamma rays.




Fig 4. Emission spectrum.




Fig 7. Image eraser lamp output.




Fig 5. Stimulated fluorescence decay. The post stimulation decay has a time constant less than 1 μ s. This shows a fraction of stored information released by red laser stimulation.




Fig 8. Effects of covering radioactive samples.

Coverings

To prevent contamination, radioactive samples are covered prior to their exposure to the Storage Phosphor Screens. There are a variety of materials used for this purpose (Fig 8).

Mylar™ provides the best transmittance for ¹⁴C particles followed by Saran™ wrap and then a hybridization bag. The thinner Mylar allows more transmittance. However, as the thickness of the Mylar decreases, the chances of screen contamination increase. Therefore, unless it is pertinent to get

as much signal as possible, a Mylar thickness of 6 µm should be used to ensure that the screen does not get contaminated. Also, the 6-µm Mylar is easier to handle than the 1.5-µm and 3-µm versions.

Physical properties

Test conditions for storage stability: 25–49 °C with 20–90% relative humidity.

Results: Storage Phosphor Screens are highly resistant to the effects of high relative humidity. Representative screens have been incubated at 32 °C/90% relative humidity for up to one year with no detectable loss in sensitivity or change in physical properties.

Test conditions for scratch resistance: 0.08 mm stylus and 25–200-g load

Results: Visible damage to the overcoat occurs when a load exceeding 50 g is applied to the stylus. Damage to the phosphor layer, seen radiographically, occurs when the load applied to the stylus exceeds 100 g.

Flatness

Screens are manufactured to be flat under conditions of 21 °C and 50% relative humidity. Actual flatness will vary with temperature, relative humidity, and screen size. The aluminum mounting plate is flat to 0.005 inches.

Asia Pacific Tel: +852 2811 8693 Fax: +852 2811 5251 **Australasia** Tel: +61 2 9899 0999 Fax: +61 2 9899 7511 **Austria** Tel: 01 576 0616 19 Fax: 01 576 0616 27
Belgium Tel: 0800 73 888 Fax: 03 272 1637 **Canada** Tel: +1 800 463 5800 Fax: +1 800 567 1008 **Central, East, and Southeast Europe** Tel: +43 1 982 3826
 Fax: +43 1 985 8327 **Denmark** Tel: 45 16 2400 Fax: 45 16 2424 **Finland & Baltics** Tel: +358 (0)9 512 39 40 Fax: +358 (0)9 512 39 439 **France** Tel: 01 69 35 67 00
 Fax: 01 69 41 96 77 **Germany** Tel: 0761 4903 490 Fax: 0761 4903 405 **Italy** Tel: 02 27322 1 Fax: 02 27302 212 **Japan** Tel: +81 3 5331 9336 Fax: +81 3 5331 9370
Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 **Middle East and Africa** Tel: +30 210 96 00 687 Fax: +30 210 96 00 693 **Netherlands** Tel: 0165 580 410
 Fax: 0165 580 401 **Norway** Tel: 2318 5800 Fax: 2318 6800 **Portugal** Tel: 21 417 70 35 Fax: 21 417 31 84 **Russia & other C.I.S. & N.I.S.** Tel: +7095 232 0250, 956 1137
 Fax: +7095 230 6377 **Southeast Asia** Tel: +60 3 8024 2080 Fax: +60 3 8024 2090 **Spain** Tel: 93 594 49 50 Fax: 93 594 49 55 **Sweden** Tel: 018 612 1900
 Fax: 018 612 1910 **Switzerland** Tel: 0848 8028 12 Fax: 0848 8028 13 **UK** Tel: 0800 616928 Fax: 0800 616927 **USA** Tel: +1 800 526 3593 Fax: +1 877 295 8102

Amersham Biosciences UK Limited

Amersham Place, Little Chalfont, Buckinghamshire, England HP7 9NA

Amersham Biosciences AB

SE-751 84 Uppsala, Sweden

Amersham Biosciences Corp

800 Centennial Avenue, PO Box 1327, Piscataway, NJ 08855 USA

Amersham Biosciences GmbH

Munzinger Strasse 9, D-79111 Freiburg, Germany

Amersham Biosciences (SV) Corp

928 East Arques Avenue, Sunnyvale CA 94085 USA

PhosphorImager, Storm, and Typhoon are trademarks of Amersham Biosciences Ltd

Amersham and Amersham Biosciences are trademarks of Amersham plc

Mylar is a trademark of DuPont.

Saran is a trademark of Dow Corporation

© Amersham Biosciences Corp. 2004 — All rights reserved

All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Biosciences group which supplies them. A copy of these terms and conditions is available on request.

www.amershambiosciences.com

