

Primer Extension Primers

Dilute the 1 µg/50 µl stocks to the following dilution for primer extension reactions:

16S rRNA Primers

23S rRNA Primers

1508	1:5	2887	1:100	1098	1:200
1490	1:100	2769*	1:100	1001	1:100
1391	1:20	2728	1:50	872	1:80
1257	1:100	2701*	1:100	806*	1:100
1199	1:50	2639	1:100	670	1:100
1046	1:60	2592*	1:50	454	1:50
906	1:50	2493	1:80	245	1:50
837	1:5	2437*	1:50	235	1:50
795	1:30	2300	1:100	<u>5S rRNA</u>	1:100
683	1:100	2274	1:100	*Nonstandard primers	
565	1:300	2145*	1:100		
480	1:30	2042	1:80		
377	1:150	1988*	1:25		
323	1:30	1905	1:50		
232	1:20	1834	1:20		
161	1:50	1679*	1:100		
		1599	1:20		
		1497*	1:100		
		1347	1:100		
		1156	1:100		